

## Preparation of mammal collected this A.M.

640314-2	<i>Microtus ochrogaster</i>	157-38-20-11-40gms ♀ ut normal
640314-3	" "	146-31-19-10-37gms ♀ ut. normal
640314-4	<i>Peromyscus maniculatus</i>	138-54-19-13-14gms ♀ ut normal
640314-5	<i>Microtus ochrogaster</i>	148-34-21-11-34gms ♀ 1x2 emb. 6mm
640314-6	" "	132-28-19.5-10-25gms ♀ ut normal
640314-7	<i>Sigmodon hispidus</i>	225-91-32-17-85gms ♂ testes 13mm
640314-8	<i>Microtus ochrogaster</i>	166-40-21-11-52gms ♀ ut normal
640314-11	" "	165-35-20-11-53gms ♂ testes 11mm
640314-13	" "	[136]-[28]-20-10-38gms ♂ testes 12mm
640314-16	<i>Sigmodon hispidus</i>	236-98-30-17-70gms ♀ ut normal
640314-21	<i>Microtus ochrogaster</i>	139-33-20-10-30gms ♂ testes 4mm
640314-23	<i>Sigmodon hispidus</i>	210-87-30-15-55gms ♀ ut. normal
640314-25	" "	295-118-34-20-130gms ♀ ut normal
640314-28	<i>Microtus ochrogaster</i>	147-33-19-10-38gms ♀ ut normal
640314-32	" "	157-36-20-10-48gms ♀ ut normal
640314-34	" "	102-24-16-8-12gms ♂ 2x2 emb. 9mm
640314-35	<i>Sigmodon hispidus</i>	218-89-31-15-64gms ♂ testes 10mm
S. only 640314-37	<i>Microtus ochrogaster</i>	162-35-20-11-45gms ♀ 2x0 emb. 12mm
640314-40	<i>Sigmodon hispidus</i>	212-82-29-15-57gms ♂ testes 10mm
640314-41	<i>Microtus ochrogaster</i>	147-33-20-10-43gms ♀ ut normal
640314-45	" "	143-33-19-11-38gms ♀ ut normal
640314-46	<i>Sigmodon hispidus</i>	245-100-32-17-82gms ♀ ut normal
640314-49	<i>Microtus ochrogaster</i>	152-33-19.5-10-40gms ♀ ut normal
640314-51	<i>Sigmodon hispidus</i>	235-98-31-16-70gms ♀ ut normal
640314-53	" "	221-87-30-18-70gms ♀ ut normal
640314-54	<i>Microtus ochrogaster</i>	148-33-20-10-40gms ♂ testes 13mm
640314-55	<i>Sigmodon hispidus</i>	225-92-31-16-65gms ♀ ut. normal
640314-57	" "	229-93-32-17-91gms ♂ testes 8mm
640314-58	<i>Microtus ochrogaster</i>	139-32-19.5-10-30gms ♂ testes 6mm
640314-62	" "	149-34-20.5-11-46gms ♀ ut. normal
640314-64	" "	156-35-20-11-50gms 1x1 emb 22mm
640314-67	" "	142-35-19-10-29gms ♂ 3x1 emb. 5mm
640314-68	<i>Sigmodon hispidus</i>	235-91-31-16-77gms ♂ testes 10mm
640314-71	<i>Microtus ochrogaster</i>	151-36-20-10-40gms ♂ testes 10mm
S. only 640314-75	" "	146-33-20-10-33gms ♂, testes 11mm
640314-80	<i>Sigmodon hispidus</i>	182-71-28-16-50gms ♀ ut normal
640314-85	" "	228-91-32-17-84gms ♂ testes 14mm
640314-88	<i>Peromyscus leucopus</i>	185-85-22.5-16-33gms ♀ ut. normal
S.K. only 640314-89	<i>Microtus ochrogaster</i>	142-31-18.5-9-25gms ♀ 3x1 emb. 18mm
640314-93	" "	153-36-20-10-42gms ♂ testes 14mm
640314-96	" "	147-31-19-10-40gms ♂ testes 12mm
640314-99	<i>Sigmodon hispidus</i>	[173]-[57]-28-16-64gms. ♀ ut normal